

Meo:PTA

PTA Catheter


Designed to Cross
in a Safer Way



Rounded tip



Tri-fold balloon wrap



OTW dual-lumen design
0.035" guide wire compatible

meo
medical

Doing Things Better

Meo:PTA

PTA Catheter

Technical Data

Catheter Design:	Over-the-wire
Shaft size:	4.8F for 5.0 – 8.0 5.8F for 9.0 – 10.0
Shaft material:	Pebax
Balloon material:	Modified PA
Balloon compliance:	Semi-compliant
Balloon fold:	Tri-fold
Cone angle:	30 degrees
Catheter tip profile:	0.067" 1.7 mm for 5.0 – 8.0 0.075" 1.9 mm for 9.0 – 10.0
Tip design:	rounded tip
Tip length:	6 mm
Nominal Pressure:	6 bar
Rated Burst Pressure:	13 (10) bar
Maximum Guide Wire:	0.035"
Usable Length / GW Lumen:	75 cm
Balloon length:	20-40-60-80 mm
Balloon diameter:	5.0-6.0-7.0-8.0-9.0-10.0 mm

Catalog Numbers

Balloon Diameter (mm)	Balloon Length (mm)			
	20	40	60	80
5.0	1-P01-5020-0	1-P01-5040-0	1-P01-5060-0	1-P01-5080-0
6.0	1-P01-6020-0	1-P01-6040-0	1-P01-6060-0	1-P01-6080-0
7.0	1-P01-7020-0	1-P01-7040-0	1-P01-7060-0	1-P01-7080-0
8.0	1-P01-8020-0	1-P01-8040-0	1-P01-8060-0	1-P01-8080-0
9.0	1-P01-9020-0	1-P01-9040-0	1-P01-9060-0	1-P01-9080-0
10.0	1-P01-10020-0	1-P01-10040-0	1-P01-10060-0	1-P01-10080-0

MeoMedical GmbH
 Potsdamer Strasse 12
 D-80802 München
 Phone: +49 (0)89-32389960
 Fax: +49 (0)89-32389970
www.meomedical.de
info@meomedical.de

distributed by

meo
medical

Doing Things Better